



Patient information: pharmacy

## Compliance aids step by step guides-Opticare®

The Opticare® is for use with round eye drop bottles. It is designed to help with the stable positioning of the bottle over your eye, so you can accurately squeeze a drop into your eye.

## Instructions for use:

- 1. After washing your hands, flip up the eye piece and open the aid. Place the neck of the bottle into the holder and close the aid.
- 2. Lie down or sit down comfortably, remove the bottle cap and clip the eye piece back down.
- 3. Tilt your head back, lower your bottom eye lid, and position the eye piece over your eye.
- 4. When ready, gently squeeze the top and bottom grooves of the aid together to release a drop.
- 5. Unless your clinician has advised otherwise, close your eye and gently press the inside corner for 30 seconds to help prevent the drop entering your mouth.
- 6. When finished, flip up the eye piece and replace the cap of your eye drop bottle.

Author: Sarah Thomas, lead satellite pharmacist and Fiona Chiu, associate chief pharmacist

Revision number: 1 Approval date: April 2018

This information can be made available in alternative formats, such as easy read or large print on request. Please call PALS: 020 7566 2324 / 020 7566 2325

