



Patient information: pharmacy

Compliance aids step by step guides- ComplEye®

The ComplEye® is specifically created for the Hylo range of multi-dose eye drops, such as Hylo-Tears and Hylo-Forte. It is designed to help you with the stable positioning of the bottle over your eye, so you can accurately squeeze a drop into your eye.

Instructions for use:

- 1. After washing your hands, place the bottle into the holder with the nozzle facing towards the eye piece.
- 2. Lie down or sit down comfortably, remove the bottle cap and tilt your head back.
- 3. With the viewing hole of the aid facing upwards, use the aid to help you lower your bottom eye lid so the lower part of the aid rests near/on your cheek bone. Place the top part of the aid on near/on your brow bone.
- 4. Position your thumbs either side of the eye piece to reach the button press at the end of the bottle.
- 5. When ready, look up through the viewing hole, and squeeze the button to release a drop.
- 6. Unless your clinician has advised otherwise, close your eye and gently press the inside corner for 30 seconds to help prevent the drop entering your mouth.

Author: Sarah Thomas, lead satellite pharmacist and Fiona Chiu, associate chief pharmacist

Revision number: 1 Approval date: April 2018

This information can be made available in alternative formats, such as easy read or large print on request. Please call PALS: 020 7566 2324 / 020 7566 2325

Moorfields Eye Hospital NHS Foundation Trust

City Road, London EC1V 2PD

Phone: 020 7253 3411 www.moorfields.nhs.uk

