



Patient information: pharmacy

Compliance aids step by step guides-Thea Eyot®

The Thea Eyot® is for use with single dose unit eye drops. It is designed to help with the stable positioning of the unit over your eye, so you can accurately squeeze a drop into your eye.

Instructions for use:

- After washing your hands, twist off the cap of the single dose unit, and gently insert it between the grooves towards the eye piece until fully in place.
- 2. Lie down or sit down comfortably and tilt your head back.
- With the viewing hole of the aid facing upwards, use the aid to help you lower your bottom eyelid, so the lower part of the aid rests near/on your cheek bone.



- 4. Place the top part of the aid on near/on your brow bone. When ready, look up through the viewing hole, and squeeze the button to release a drop.
- 5. Unless your clinician has advised otherwise, close your eye and gently press the inside corner for 30 seconds to help prevent the drop entering your mouth.

Author: Sarah Thomas, lead satellite pharmacist and Fiona Chiu, associate chief pharmacist

Revision number: 1 Approval date: April 2018

This information can be made available in alternative formats, such as easy read or large print on request. Please call PALS: 020 7566 2324 / 020 7566 2325